

CardioWare™ Harness

Diagnostic quality ECG without skin preparation or exposed wires

SILVER™ bumps



Orbital Research Inc.

Creating tomorrow's products from today's innovations.

*"The **Bumps** have it when it comes to CardioWare." Orbital's SilverBumps electrodes are key to 'locking' sensors to the skin.*

Medical Design Magazine, June 6, 2013



Key Features:


- Equipped with SilverBumps™
- Adjusts to fit most body types
- No exposed wires
- Disposable

Benefits:

- Eliminate exposed wires
- Minimize skin preparation
- Don/doff quickly and easily
- Improve patient comfort
- Increase patient compliance
- Extend monitoring time

Health Applications:

- Cardiac Monitoring
- Mobile Outpatient Telemetry
- Disease Management
- Preventative Medicine
- Elite Sport/Athlete Training



Two ECG plots attained using the CardioWare Harness from a test subject with a known arrhythmia. The CardioWare harness was able to clearly identify the arrhythmia.

Orbital Research, Inc.

4415 Euclid Avenue, Suite 500
Cleveland, Ohio 44103-3757

Phone: 216-649-0399

Fax: 216-649-0347

E-mail: Info@orbitalresearch.com

www.silverbumps.com